

# DERMA PEEL®

[www.genesisbiosystems.com](http://www.genesisbiosystems.com)



THE LEADER IN MICRODERMABRASION

Visit our website  
for more information, videos, and testimonials  
about our other innovative products.

Our goal is to provide you with proven tools that  
will assist you in your success.

**Genesis**  
Biosystems  
INNOVATIONS IN REGENERATION

© 2012 Genesis Biosystems, Inc.  
DRP-BRO-001 Rev A

Genesis Biosystems, Inc.  
1500 Eagle Court  
Lewisville, TX 75057  
(888) 577-PEEL (7335)  
(972) 315-7888

# DERMA PEEL®

Unlike any other.

State-of-the-art  
**crystal-free**  
**microdermabrasion**  
with proven clinical results, high patient  
satisfaction and zero maintenance.

FROM THE LEADER IN  
MICRODERMABRASION

**Genesis**  
Biosystems  
INNOVATIONS IN REGENERATION





State-of-the-art crystal-free microdermabrasion with proven clinical results, high patient satisfaction and zero maintenance.



### It's a fact.

More people request microdermabrasion than any other aesthetic procedure.

Trust the brand that more than 5,000 plastic surgery and dermatology clinics throughout the world have chosen - **Genesis Biosystems.**

### Why Choose Derma Peel®?

#### Crystal-Free Technology

Crystal-free microdermabrasion gives a superior patient and practitioner experience with its safe, even and reproducible results.

#### Saves You Time

Derma Peel® utilizes disposable crystal-free tips so there's no need to scrub, soak or sterilize - simply discard and move on to your next client.

#### Safety First

Our no-maintenance system provides a safer experience for you and the patient. There is no used crystal to dispose and no hazardous biomedical waste.

#### Dependability

Derma Peel® is built to last with a powerful, efficient motor that will sustain you through years of successful treatments and ongoing revenue.

#### We've Got You Covered

Derma Peel's® three year warranty covers all parts and labor. Should your machine need service, we provide a loaner unit so you stay up and running.

#### Training and Certification

We offer thorough, online training and certification so that the entire staff feels confident and comfortable operating the system.

### How does Derma Peel® work?

Derma Peel® utilizes crystal-free tips that rely on constant vacuum pressure to perform an even and effective peel, rather than requiring the practitioner to apply pressure. This makes Derma Peel® a much more comfortable procedure for the practitioner to provide and gives your patient outstanding results.



Our Bio-Barrier™ disposable tip design eliminates concerns about cross-contamination. The abraded surface is completely surrounded by vacuum, ensuring that all cells and bacteria are whisked away, leaving the skin clean and healthy.



Three grit strengths (240, 180 and 120) give you the flexibility to offer various depths of exfoliation.

### Microdermabrasion has proven to be clinically effective for:

- Dry or Oily Skin • Fine Lines & Wrinkles
- Sun Damaged Skin • Hyperpigmentation • Acne Prone Skin
- Comedone Removal • Scar Tissue Reduction



Microdermabrasion is a non-invasive, painless procedure that treats all skin types. Patients experience zero downtime, immediate results and instant gratification.

### Why is Genesis Biosystems the best-selling microdermabrasion brand?



- State-of-the-art technology
- Completely maintenance free
- No cleaning or sterilizing
- Excellent, reproducible results
- Recurring patient satisfaction